


# VALUE-LINE Ductless Fume Hoods

900mm or 1200mm Units



## VALUE-LINE Ductless Fume Hoods, 900mm or 1200mm Units

The VALUE-LINE Ductless Fume Hoods, constructed from fire-retardant PPS (polyphenylene sulfide) material, ensure a safe working environment for laboratory staff handling various chemicals.


The structure of the fume cabinet includes 8mm welded PPS, which offers high chemical resistance. The hood also features a 5mm safety front glass window with a sliding front sash counterbalanced by a counterweight, and side windows to improve visibility.

The fume hood is equipped with a 7-inch colour touchscreen display, and an airflow alarm and display. A red chamber alert light signals critical alarms. Additional features include an environmentally friendly EC fan with 10 levels of speed adjustment, 800 lux LED lighting, and adjustable air velocity at  $100 \pm 20$  fpm ( $0.5 \pm 0.1$  m/s).

The fume hood offers back access for filter replacement and includes a 250mm exhaust point collar for connection to an external exhaust system. It is certified to meet CE standards.



Carbon Filter



### Product Features

- Adjustable Air Velocity of  $100 \pm 20$  fpm ( $0.5 \pm 0.1$  m/s)
- Tempered Glass Side Windows for Improved Visibility
- Red Chamber Alert Alarm for Critical Alarms
- 10-Speed Environmentally Friendly EC Fan
- Back Access for Filter Replacement

### Specifications

#### Product Information

Type	VALUE-LINE Ductless Fume Hoods
Hood Material	PPS Structure & Built-In Sealed PPS Worktop
Digital Display	7" Colour Touchscreen Display
Illumination	800 lux, Eco-Friendly LED Lighting
Microprocessor Control	7" Colour Touchscreen, Airflow System, 10-Speed Fan Control, Filter Replacement Alarm, Service Alarm & Lights Control

**Publication Date:** 21/08/2024

**Westlab Pty. Ltd.**  
4 Cargo Way, Mitchell Park,  
Vic 3355, Australia

**Australia**  
1800 358 101  
sales.au@westlab.com  
www.westlab.com.au

**Canada**  
1 877 822 1455  
sales.ca@westlab.com  
www.westlab.com

**New Zealand**  
+64 9 553 3677  
sales.nz@westlab.com  
www.westlab.co.nz

**United States**  
646 568 5391  
sales.us@westlab.com  
www.westlab.com

### Performance Information

Adjustable Air Velocity	0.5 ± 0.1m/s, 100 ± 20 fpm
Power Supply	110/220 V, 50/60 Hz, Single Phase
Noise Level	< 55 dB (900mm Unit) or < 57 dB (1200mm Unit) <i>Tested 200mm from the workstation and 1.2M above the ground.</i>
EC Fan	Yes, 10-Speed Fan Control

### Physical Information

Front Window	5mm safety front glass window with a sliding front sash with a counterweight.
Front Sash Opening	500mm
Side Windows	Yes, Tempered Glass
Exhaust Point Collar	250mm collar for connection to an external exhaust.

### Safety Information

Chamber Alert Light	A red light for critical alarms.
Country of Origin	United States of America
Standards Conformity	CE

### Product Information – 900mm Unit

Code	665-601
Name	VALUE-LINE 900 Ductless Fume Hood, PPS
Workspace Dimensions	880mm W x 450mm D x 550mm H
External Dimensions	900mm W x 630mm D x 1050mm H
Net Weight	100 kg

### Product Information – 1200mm Unit

Code	665-602
Name	VALUE-LINE 1200 Ductless Fume Hood, PPS
Workspace Dimensions	1180mm W x 450mm D x 550mm H
External Dimensions	1200mm W x 630mm D x 1050mm H
Net Weight	120 kg

**Australia**  
1800 358 101  
sales.au@westlab.com  
www.westlab.com.au


**Canada**  
1 877 822 1455  
sales.ca@westlab.com  
www.westlab.com


**New Zealand**  
+64 9 553 3677  
sales.nz@westlab.com  
www.westlab.co.nz

**United States**  
646 568 5391  
sales.us@westlab.com  
www.westlab.com


## Bundled Configurations


### 900mm Bundles

<b>Code</b>	<b>665-611</b>	
<b>Description</b>	VALUE-LINE 900 Ductless Fume Hood, Carbon Filter & Metal Stand	
<b>This bundle includes</b>		
<b>Title</b>	<b>Qty</b>	
VALUE-LINE 900 Ductless Fume Hood, PPS, 900 W x 630 D x 1050mm H	1	
Carbon Filter for 900 VALUE-LINE Ductless Fume Hood	1	
Metal Stand & Castors for 900 VALUE-LINE Ductless Fume Hood	1	

<b>Code</b>	<b>665-612</b>	
<b>Description</b>	VALUE-LINE 900 Ductless Fume Hood & Carbon Filter	
<b>This bundle includes</b>		
<b>Title</b>	<b>Qty</b>	
VALUE-LINE 900 Ductless Fume Hood, PPS, 900 W x 630 D x 1050mm H	1	
Carbon Filter for 900 VALUE-LINE Ductless Fume Hood	1	

### 1200mm Bundles

<b>Code</b>	<b>665-613</b>	
<b>Description</b>	VALUE-LINE 1200 Ductless Fume Hood, Carbon Filter & Metal Stand	
<b>This bundle includes</b>		
<b>Title</b>	<b>Qty</b>	
VALUE-LINE 1200 Ductless Fume Hood, PPS, 1200 W x 630 D x 1050mm H	1	
Carbon Filter for 1200 VALUE-LINE Ductless Fume Hood	1	
Metal Stand & Castors for 1200 VALUE-LINE Ductless Fume Hood	1	

<b>Code</b>	<b>665-614</b>	
<b>Description</b>	VALUE-LINE 1200 Ductless Fume Hood & Carbon Filter	
<b>This bundle includes</b>		
<b>Title</b>	<b>Qty</b>	
VALUE-LINE 1200 Ductless Fume Hood, PPS, 1200 W x 630 D x 1050mm H	1	
Carbon Filter for 1200 VALUE-LINE Ductless Fume Hood	1	

**Australia**  
1800 358 101  
sales.au@westlab.com  
www.westlab.com.au

**Canada**  
1 877 822 1455  
sales.ca@westlab.com  
www.westlab.com

**New Zealand**  
+64 9 553 3677  
sales.nz@westlab.com  
www.westlab.co.nz

**United States**  
646 568 5391  
sales.us@westlab.com  
www.westlab.com

## Ordering Information

SKU	Description
665-601	VALUE-LINE 900 Ductless Fume Hood, PPS, 900 W x 630 D x 1050mm H
665-602	VALUE-LINE 1200 Ductless Fume Hood, PPS, 1200 W x 630 D x 1050mm H
665-604	Metal Stand & Castors for 900 VALUE-LINE Ductless Fume Hood
665-605	Metal Stand & Castors for 1200 VALUE-LINE Ductless Fume Hood
665-607	Carbon Filter for 900 VALUE-LINE Ductless Fume Hood
665-608	Carbon Filter for 1200 VALUE-LINE Ductless Fume Hood



### Australia

1800 358 101  
[sales.au@westlab.com](mailto:sales.au@westlab.com)  
[www.westlab.com.au](http://www.westlab.com.au)

### Canada

1 877 822 1455  
[sales.ca@westlab.com](mailto:sales.ca@westlab.com)  
[www.westlab.com](http://www.westlab.com)

### New Zealand

+64 9 553 3677  
[sales.nz@westlab.com](mailto:sales.nz@westlab.com)  
[www.westlab.co.nz](http://www.westlab.co.nz)

### United States

646 568 5391  
[sales.us@westlab.com](mailto:sales.us@westlab.com)  
[www.westlab.com](http://www.westlab.com)